

# Rat VEGFR2/Flk-1/CD309/KDR (His Tag) recombinant protein



Catalog Number: 504258

## General Information

### Protein Construction

A DNA sequence encoding the rat KDR (AAH87029.1) extracellular domain (Met 1-Glu 760) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Rat

### Expression Host

Human Cells

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Ala 20

### Molecular Mass

The recombinant rat KDR comprises 752 amino

acids and predicts a molecular mass of 84.3 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Rat VEGFR2 / Flk-1 / CD309 / KDR Protein (His Tag) SDS-PAGE