

# Canine CD40/TNFRSF5 (Fc Tag) recombinant protein



Catalog Number: 504279

## General Information

### Gene Name Synonym

B-cell surface antigen CD40; CD40L receptor

### Protein Construction

A DNA sequence encoding the canine CD40 (Q7YRL5) (Met1-Ala194) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Canine

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Glu 21

### Molecular Mass

The recombinant canine CD40/Fc is a disulfide-linked homodimer. The reduced monomer

comprises 415 amino acids and has a predicted molecular mass of 46.1 kDa. The apparent molecular mass of the protein is approximately 50-56 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Canine CD40 / TNFRSF5 Protein (Fc Tag) SDS-PAGE