

# H5N1 HA (His Tag) recombinant protein



Catalog Number: 504384

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/turkey/Turkey/1/2005(H5N1)) hemagglutinin (Q207Z6, with aa 339-342 mutated from "GERR" to "RESR") (Met1-Gln531), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H5N1

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 17

### Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/turkey/Turkey/1/2005(H5N1)) consists

526 amino acids and predicts a molecular mass of 60.1 kDa.

### Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 %glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA Protein (His Tag)