

# H10N8 HA (full length) recombinant protein



Catalog Number: 504466

## General Information

### Protein Construction

A DNA sequence encoding the full length of Influenza A virus (A/Jiangxi-Donghu/346/2013(H10N8)) hemagglutinin (translated amino acid sequence of EPI497477 ) (Met1-Ile561), termed as HA, was expressed.

### Organism

H10N8

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Measured by its ability to agglutinate guinea pig red blood cells.

HA titer is 30-150 $\mu$ g/mL for 1%GRBC

### Purity

> 90 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Leu 17

### Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Jiangxi-Donghu/346/2013(H10N8)) consists 545 amino acids and predicts a molecular mass of 60.6 kDa.

### Formulation

Lyophilized from sterile 20mM Tris, 150mM NaCl, pH7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H10N8 (A/Jiangxi-Donghu/346/2013) Hemagglutinin / HA0 Protein (full length)