

# Rabbit TrkA/NTRK1 (His Tag) recombinant protein



Catalog Number: 504479

---

## General Information

### Protein Construction

A DNA sequence encoding the rabbit NTRK1 (XP\_008262512.1) (Ala33-Glu414) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Rabbit

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

## Predicted N terminal

Ala 33

## Molecular Mass

The recombinant rabbit NTRK1 consists of 393 amino acids and predicts a molecular mass of 43.4 kDa.

## Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

## Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

## Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.