

Canine KIM-1/TIM1/HACVR1 (His Tag) recombinant protein



Catalog Number: 504557

General Information

Protein Construction

A DNA sequence encoding the extracellular domain of canine HAVCR1 (XP_854888.1) (Tyr 21-Ser 160) was expressed with a polyhistidine tag at the N-terminus.

Organism

Canine

Expression Host

Human Cells

QC Testing

Purity

> 96 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant canine HAVCR1 consists of 156 amino acids and has a calculated molecular mass of 17.6 kDa. In SDS-PAGE under reducing

conditions, the apparent molecular mass of the protein is approximately 35-40 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Canine KIM-1 / TIM1 / HAVCR1 Protein (His Tag)
SDS-PAGE