

H3N2 HA (His Tag) recombinant protein



Catalog Number: 504573

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met1-Arg345) of the influenza A hemagglutinin (A/Wisconsin/67/2005 (H3N2)) (ACF54576.1), termed as HA1, was expressed, fused with a C-terminal polyhistidine tag.

Organism

H3N2

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The secreted recombinant influenza A H3N2 HA1 subunit (A/Wisconsin/67/2005 (H3N2)) comprises

340 amino acids and has a predicted molecular mass of 38 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H3N2 (A/Wisconsin/67/X-161/2005) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE