

# Canine RBP4 (His Tag) recombinant protein



Catalog Number: 504675

## General Information

### Protein Construction

A DNA sequence encoding the canine RBP4 (F1PQI4) (Met1-Leu201) was expressed with a C-terminal polyhistidine tag.

### Organism

Canine

### Expression Host

Human Cells

## QC Testing

### Activity

Measured by its ability to bind alltrans retinoic acid. The binding of retinoic acid results in the quenching of Trp fluorescence in RBP4. The 50% binding concentration (BC50) is  $>1 \mu\text{M}$

### Purity

$> 94 \%$  as determined by SDS-PAGE

### Endotoxin

$< 1.0 \text{ EU per } \mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^\circ\text{C}$

### Predicted N terminal

Glu 19

### Molecular Mass

The recombinant canine RBP4 comprises 194 amino acids and has a predicted molecular mass of 22.4 kDa. The apparent molecular mass of the protein is approximately 23 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^\circ\text{C}$  to  $-80^\circ\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Canine RBP4 Protein (His Tag) SDS-PAGE