

# Sus scrofa (Pig) Legumain (His Tag) recombinant protein



Catalog Number: 504680

## General Information

### Protein Construction

A DNA sequence encoding the sus scrofa(Pig) LGMN (XP\_001927117.4) (Met1-Tyr433) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Sus scrofa (Pig)

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Val 18

### Molecular Mass

The recombinant sus scrofa (Pig) LGMN consists of 427 amino acids and has a predicted molecular mass of 48.8 kDa. The apparent molecular mass of

it is approximately 53.7 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Sus scrofa (Pig) Legumain / LGMN / AEP Protein (His Tag) SDS-PAGE