

Rat VE-Cadherin/CD144 (Fc Tag) recombinant protein



Catalog Number: 504745

General Information

Protein Construction

A DNA sequence encoding the rat CDH5 (NP_001100877.1) (Met1-Gln585) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 75 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 25

Molecular Mass

The recombinant rat CDH5 comprises 802 amino acids and predicts a molecular mass of 90.7 kDa. The apparent molecular mass of the recombinant

protein is approximately 116 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat VE-Cadherin / CD144 / CDH5 Protein (Fc Tag) SDS-PAGE