

Mouse LAG3 / CD223 / Lymphocyte activation gene 3 (His Tag) recombinant protein



Catalog Number: 504764

General Information

Gene Name Synonym:

Protein Construction:

A DNA sequence encoding the mouse Lag3 (NP_032505.1) (Met1-Leu442) was expressed with a polyhistidine tag at the C-terminus.

Source: Mouse

Expression Host: Human Cells

QC Testing

Activity:

Purity:

> 95 % as determined by SDS-PAGE.

Endotoxin:

<1.0 EU per μg protein as determined by the LAL method.

Stability:

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal:

Ser 23

Molecular Mass:

The recombinant mouse Lag3 consists 431 amino

acids and predicts a molecular mass of 46.8 kDa.

Formulation:

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage:


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution:

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE:



Mouse LAG3 / CD223 / Lymphocyte activation gene 3 Protein (His Tag)