

Human Amyloid beta A4 protein (aa 672-710, His & GST Tag)



Catalog Number: 504849

General Information

Protein Construction

A DNA sequence encoding the amino acids (Asp672-Val710) of human Amyloid beta A4 protein (APP770) (P05067-1), corresponding to the Beta-amyloid protein 39, was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human Beta-APP39/GST chimera consists of 278 amino acids and has a calculated molecular mass of 32.2 KDa. It migrates as an

approximately 34 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50 mM Tris, 100 mM NaCl, 20% glycerol, 0.05% Tween 20, pH 9.5.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human Beta-amyloid 39 / Beta-APP39 Protein (aa 672-710, His & GST Tag) SDS-PAGE