

Canine CD122 / IL-2RB (Fc Tag) recombinant protein



Catalog Number: 504871

General Information

Protein Construction

A DNA sequence encoding the extracellular domain of the canine IL2RB (XP_003431552.1) (Met1-Pro244) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Canine

Expression Host

Human Cells

QC Testing

Purity

> 80 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 25

Molecular Mass

The recombinant canine IL2RB/Fc is a disulfide-linked homodimer. The reduced monomer comprises 461 amino acids and has a predicted

molecular mass of 52.3 kDa. The apparent molecular mass of the protein is approximately 67 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Canine IL2RB / IL2 Receptor beta Protein (Fc Tag) SDS-PAGE