

Cynomolgus EGFR/HER1/ErbB1 (His Tag) recombinant protein



Catalog Number: 504912

General Information

Protein Construction

A DNA sequence encoding the cynomolgus EGFR (XP_005549616.1) (Met1-Ser645) was expressed with a polyhistidine tag at the C-terminus.

Organism

Cynomolgus

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 25

Molecular Mass

The recombinant cynomolgus EGFR consists of 632 amino acids and predicts a molecular mass of 70.1 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus EGFR / HER1 / ErbB1 Protein (His Tag)