

General Information

Gene Name Synonym

Thy-1 membrane glycoprotein variant 2; cDNA FLJ33325 fis, clone BNGH42008832, highly similar to THY-1 MEMBRANE GLYCOPROTEIN; cDNA FLJ33677 fis, clone BRAWH2002125, highly similar to THY-1 MEMBRANE GLYCOPROTEIN

Protein Construction

A DNA sequence encoding the human THY1 (NP_006279.2) (Met1-Cys130) was expressed with five amino acids (DDDDK) at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 20

Molecular Mass

The recombinant human THY1 consists of 117 amino acids and predicts a molecular mass of 13.2 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

SDS-PAGE



Human CD90/THY-1 Protein 15092