

# Rat SerpinG1 recombinant protein (C-His)



Catalog Number: 514813

## General Information

### Gene Name Synonym

C1 Inh; C1Inh; C1 esterase inhibitor; C1-inhibiting factor; Serpin G1

### Protein Construction

A DNA sequence encoding the rat SerpinG1 (NP\_954524.1) (Met1-Ala504) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Rat

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 23

### Molecular Mass

The recombinant rat SerpinG1 consists of 493 amino acids and predicts a molecular mass of 54.7 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.


### Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

### SDS-PAGE



Rat C1 inhibitor/SerpinG1 Protein 15531



Rat C1 inhibitor / SerpinG1 Protein (His Tag)