

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Philippines/472/2002 (MDCK)(H3N2)) hemagglutinin (AIU46063.1) (Met1-Trp530) was expressed with a polyhistidine tag at the C-terminus.

Organism

H3N2

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95% as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The recombinant HA subunit of the Influenza A virus (A/Philippines/472/2002 (MDCK)(H3N2))

hemagglutinin consists of 525 amino acids and predicts a molecular mass of 59.1 kDa.

Formulation

Lyophilized from sterile 20mM PBS Ph7.0, 500mM NaCl, 10% gly.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

SDS-PAGE



H3N2 HA Protein 15911