

DENV DENV-E recombinant protein (C-His)



Catalog Number: 514893

General Information

Protein Construction

A DNA sequence encoding the envelope protein of Dengue virus (DENV)(type 4, strain Philippines/H241/1956) (AAX48017.1) (Met280-Gly674) was expressed with a polyhistidine tag at the C-terminus.

Organism

DENV

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met 280

Molecular Mass

The recombinant envelope protein of Dengue virus (DENV)(type 4, strain Philippines/H241/1956) consists of 406 amino acids and predicts a

molecular mass of 44.8 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 %gly, pH 8.0.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

SDS-PAGE



Dengue virus (DENV) (type 4, strain Philippines/H241/1956) E / Envelope Protein (ECD, His Tag)