

Human GLCE/KIAA0836 (Flag & His tag) recombinant protein



Catalog Number: 518029

General Information

Gene Name Synonym

D-glucuronyl C5-epimerase; Heparan sulfate C5-epimerase; Hsepi; Heparin/heparan sulfate:glucuronic acid C5-epimerase; Heparosan-N-sulfate-glucuronate 5-epimerase; GLCE

Protein Construction

A DNA sequence encoding the Human GLCE/KIAA0836 Asn29-Asn617 is expressed with a 6His, Flag tag at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20 mM Tris, 200mM NaCl, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human GLCE/KIAA0836 (Flag & His tag)
recombinant protein