

Human PRTN3/MBN (His tag) recombinant protein



Catalog Number: 518034

General Information

Gene Name Synonym

Myeloblastin, AGP7, C-ANCA Antigen, Leukocyte Proteinase 3, PR-3, PR3, Neutrophil Proteinase 4, NP-4, P29, Wegener Autoantigen, PRTN3, MBN

Protein Construction

A DNA sequence encoding the Human PRTN3/MBN Ala26-Arg249 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 10mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PRTN3/MBN (His tag) recombinant protein