

# Human PSMD10 recombinant protein



Catalog Number: 518111

## General Information

### Gene Name Synonym

26S Proteasome Non-ATPase Regulatory Subunit 10, 26S Proteasome Regulatory Subunit p28, Gankyrin, p28(GANK), PSMD10

### Protein Construction

A DNA sequence encoding the Human PSMD10 Met1-Gly226 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PSMD10 recombinant protein