

Mouse Sema4d/Semacl2/Semaj (Fc tag) recombinant protein



Catalog Number: 518149

General Information

Gene Name Synonym

SEMA4D/Semaphorin-4D/M-Sema G/Semaphorin-C-like 2/Semaphorin-J/Sema J/CD100/Semacl2/Semaj

Protein Construction

A DNA sequence encoding the Mouse Sema4d/Semacl2/Semaj Phe24-Met711 is expressed with a Fc tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Mouse Sema4d/Semacl2/Semaj (Fc tag)
recombinant protein