

Human SULT1A2/STP2 (His tag) recombinant protein



Catalog Number: 518157

General Information

Gene Name Synonym

Sulfotransferase 1A2, ST1A2, Aryl Sulfotransferase 2, Phenol Sulfotransferase 2, Phenol-Sulfating Phenol Sulfotransferase 2, P-PST 2, SULT1A2, STP2

Protein Construction

A DNA sequence encoding the Human SULT1A2/STP2 Met1-Leu295 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human SULT1A2/STP2 (His tag) recombinant protein