

# Human SULT1A2/STP2 (His tag) recombinant protein



Catalog Number: 518157

## General Information

### Gene Name Synonym

Sulfotransferase 1A2, ST1A2, Aryl  
Sulfotransferase 2, Phenol Sulfotransferase 2,  
Phenol-Sulfating Phenol Sulfotransferase 2, P-PST  
2, SULT1A2, STP2

### Protein Construction

A DNA sequence encoding the Human  
SULT1A2/STP2 Met1-Leu295 is expressed with a  
6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM  
Tris, 100mM NaCl, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SULT1A2/STP2 (His tag) recombinant  
protein