

Human PDIA4/ERP70/ERP72 (His tag) recombinant protein



Catalog Number: 518164

General Information

Gene Name Synonym

Protein Disulfide-Isomerase A4, Endoplasmic Reticulum Resident Protein 70, ER Protein 70, ERp70, Endoplasmic Reticulum Resident Protein 72, ER Protein 72, ERp-72, ERp72, PDIA4, ERP70, ERP72

Protein Construction

A DNA sequence encoding the Human PDIA4/ERP70/ERP72 Val21-Leu645 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 150mM NaCl, 10%Glycerol, pH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PDIA4/ERP70/ERP72 (His tag)
recombinant protein