

# Human PDIA4/ERP70/ERP72 (His tag) recombinant protein



Catalog Number: 518164

## General Information

### Gene Name Synonym

Protein Disulfide-Isomerase A4, Endoplasmic Reticulum Resident Protein 70, ER Protein 70, ERp70, Endoplasmic Reticulum Resident Protein 72, ER Protein 72, ERp-72, ERp72, PDIA4, ERP70, ERP72

### Protein Construction

A DNA sequence encoding the Human PDIA4/ERP70/ERP72 Val21-Leu645 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PDIA4/ERP70/ERP72 (His tag)  
recombinant protein