

Human GABARAPL1/GEC1 (His tag) recombinant protein



Catalog Number: 518166

General Information

Gene Name Synonym

Gamma-Aminobutyric Acid Receptor-Associated Protein-Like 1, Early Estrogen-Regulated Protein, GABA(A) Receptor-Associated Protein-Like 1, Glandular Epithelial Cell Protein 1, GEC-1, GABARAPL1, GEC1

Protein Construction

A DNA sequence encoding the Human GABARAPL1/GEC1 Met1-Lys117 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human GABARAPL1/GEC1 (His tag) recombinant protein