

# Human S100P/S100E (His tag) recombinant protein



Catalog Number: 518185

## General Information

### Gene Name Synonym

Protein S100-P, Protein S100-E, S100 Calcium-Binding Protein P, S100P, S100E

### Protein Construction

A DNA sequence encoding the Human S100P/S100E Met1-Lys95 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM Tris, 1mM DTT, 0.1M NaCl, 20% Glycerol, pH 8.0 .

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

## SDS-PAGE



Human S100P/S100E (His tag) recombinant protein