

Human UBE2M/UBC12 recombinant protein



Catalog Number: 518186

General Information

Gene Name Synonym

NEDD8-conjugating enzyme Ubc12, NEDD8 carrier protein, NEDD8 protein ligase, Ubiquitin-conjugating enzyme E2 M, UBC12, UBE2M,

Protein Construction

A DNA sequence encoding the Human UBE2M/UBC12 Met1-Lys183 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, pH 7.5, 2mM DTT, 150mM NaCl, 10% glycerol .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2M/UBC12 recombinant protein