

Mouse Cst7 (His tag) recombinant protein



Catalog Number: 518190

General Information

Gene Name Synonym

Cystatin-F, Cystatin-like Metastasis-Associated Protein, Leukocystatin, CMAP, Cystatin-7, Cst7,

Protein Construction

A DNA sequence encoding the Mouse Cst7 Ala19-Gln144 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Cst7 (His tag) recombinant protein