

# Human IL20RB/DIRS1/UNQ557/PRO1114 (Fc tag) recombinant protein



Catalog Number: 518192

## General Information

### Gene Name Synonym

Interleukin-20 receptor subunit beta;IL-20 receptor subunit beta;IL-20R-beta;IL-20RB;IL-20R2;DIRS1;hCG\_2022374;FNDC6; MGC34923; fibronectin type III domain containing 6; interleukin-20 receptor II

### Protein Construction

A DNA sequence encoding the Human IL20RB/DIRS1/UNQ557/PRO1114 Asp30-Ala230 is expressed with a Fc tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

## Formulation

Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4 .

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human IL20RB/DIRS1/UNQ557/PRO1114 (Fc tag)  
recombinant protein