

# Human HAO1/GOX1/HAOX1 (Trx & His tag) recombinant protein



Catalog Number: 518226

## General Information

### Gene Name Synonym

Hydroxyacid Oxidase 1, HAOX1, Glycolate Oxidase, GOX, HAO1, GOX1, HAOX1

### Protein Construction

A DNA sequence encoding the Human HAO1/GOX1/HAOX1 Met1-Ile370 is expressed with a Trx, 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

## SDS-PAGE



Human HAO1/GOX1/HAOX1 (Trx & His tag)  
recombinant protein