

# Mouse Cpq/Hls2/Pgcp (His tag) recombinant protein



Catalog Number: 518228

## General Information

### Gene Name Synonym

PGCP/Carboxypeptidase Q/Hematopoietic lineage switch 2/Plasma glutamate carboxypeptidase/Cpq/Hls2

### Protein Construction

A DNA sequence encoding the Mouse Cpq/Hls2/Pgcp Lys19-Ser470 is expressed with a 6His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Cpq/Hls2/Pgcp (His tag) recombinant protein