

# Human UBE2D3/UBC5C/UBCH5C recombinant protein



Catalog Number: 518237

## General Information

### Gene Name Synonym

Ubiquitin-conjugating enzyme E2 D3, Ubiquitin carrier protein D3, Ubiquitin-conjugating enzyme E2(17)KB 3, Ubiquitin-conjugating enzyme E2-17 kDa 3, Ubiquitin-protein ligase D3, UBE2D3 and UBCH5C.

### Protein Construction

A DNA sequence encoding the Human UBE2D3/UBC5C/UBCH5C Met1-Met147 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% glycerol, pH 7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human UBE2D3/UBC5C/UBCH5C recombinant protein