

# Human RNASET2/RNASE6PL (His tag) recombinant protein



Catalog Number: 518242

## General Information

### Gene Name Synonym

Ribonuclease T2,3.1.27.-,Ribonuclease  
6,RNASE6PL

### Protein Construction

A DNA sequence encoding the Human  
RNASET2/RNASE6PL Asp25-His256 is expressed  
with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL  
method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM  
TrisHcl, 150mM NaCl,20%Glycerol,pH7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

## SDS-PAGE



Human RNASET2/RNASE6PL (His tag)  
recombinant protein