

Human RNASET2/RNASE6PL (His tag) recombinant protein



Catalog Number: 518242

General Information

Gene Name Synonym

Ribonuclease T2,3.1.27.-,Ribonuclease
6,RNASE6PL

Protein Construction

A DNA sequence encoding the Human
RNASET2/RNASE6PL Asp25-His256 is expressed
with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM
TrisHcl, 150mM NaCl,20%Glycerol,pH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human RNASET2/RNASE6PL (His tag)
recombinant protein