

Human LALBA/LYZL7 (His tag) recombinant protein



Catalog Number: 518249

General Information

Gene Name Synonym

Lactose synthase B protein, Lysozyme-like protein 7, LALBA, LYZL7

Protein Construction

A DNA sequence encoding the Human LALBA/LYZL7 Lys20-Leu142 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM Tris, 150mM NaCl, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human LALBA/LYZL7 (His tag) recombinant protein