

Human UBE2Z/HOYS7 (His tag) recombinant protein



Catalog Number: 518275

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 Z, Uba6-Specific E2 Conjugating Enzyme 1, Use1, Ubiquitin Carrier Protein Z, Ubiquitin-Protein Ligase Z, UBE2Z, HOYS7

Protein Construction

A DNA sequence encoding the Human UBE2Z/HOYS7 Met1-Val246 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 50mM NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2Z/HOYS7 (His tag) recombinant protein