

Human PDIA5 (His tag) recombinant protein



Catalog Number: 518307

General Information

Gene Name Synonym

Protein disulfide-isomerase A5, Protein disulfide isomerase-related protein, PDIA5, PDIR

Protein Construction

A DNA sequence encoding the Human PDIA5 Ser22-Leu262 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of PBS,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PDIA5 (His tag) recombinant protein