

# Human S100A8/CAGA/CFAG/MRP8 (His tag) recombinant protein



Catalog Number: 518312

## General Information

### Gene Name Synonym

Protein S100-A8, S100A8, Calgranulin-A, Cystic fibrosis antigen, Leukocyte L1 complex light chain, MRP-8

### Protein Construction

A DNA sequence encoding the Human S100A8/CAGA/CFAG/MRP8 Met1-Glu93 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

## Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of 20mM HEPES, 150mM NaCl, pH 7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human S100A8/CAGA/CFAG/MRP8 (His tag)  
recombinant protein