

Human YKT6 (His tag) recombinant protein



Catalog Number: 518319

General Information

Gene Name Synonym

Synaptobrevin homolog YKT6, YKT6

Protein Construction

A DNA sequence encoding the Human YKT6 Met1-Met198 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 50mM Tris, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human YKT6 (His tag) recombinant protein