

# Human XPNPEP1/XPNPEPL/XPNPEPL1 (His tag) recombinant protein



Catalog Number: 518334

## General Information

### Gene Name Synonym

Xaa-Pro Aminopeptidase 1, Aminoacylproline Aminopeptidase, Cytosolic Aminopeptidase P, Soluble Aminopeptidase P, sAmp, X-Pro Aminopeptidase 1, X-Prolyl Aminopeptidase 1 Soluble, XPNPEP1, XPNPEPL, XPNPEPL1

### Protein Construction

A DNA sequence encoding the Human XPNPEP1/XPNPEPL/XPNPEPL1 Pro2-His623 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 10% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human XPNPEP1/XPNPEPL/XPNPEPL1 (His tag)  
recombinant protein