

Mouse Il1rl1/Ly84/St2/Ste2 (Fc tag) recombinant protein



Catalog Number: 518342

General Information

Gene Name Synonym

Interleukin-1 receptor-like 1, Lymphocyte antigen 84, Protein ST2, Protein T1, Il1rl1, DER4, Fit-1, IL-1 R4, IL-1RL1, IL-33R, Ly84, ST2L, ST2V, suppression of tumorigenicity 2

Protein Construction

A DNA sequence encoding the Mouse Il1rl1/Ly84/St2/Ste2 Ser27-Ala337 is expressed with a Fc tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS, pH 7.4.

Usage Guide


Storage

Lyophilized protein should be stored at $<-20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Mouse Il1rl1/Ly84/St2/Ste2 (Fc tag) recombinant protein