

Human BCL2L1/BCL2L/BCLX (His tag) recombinant protein



Catalog Number: 518346

General Information

Gene Name Synonym

Bcl-2-Like Protein 1, Bcl2-L-1, Apoptosis Regulator
Bcl-X, BCL2L1, BCL2L, BCLX

Protein Construction

A DNA sequence encoding the Human
BCL2L1/BCL2L/BCLX Met1-Arg212 is expressed
with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM
HEPES, 50mM KCl, 20% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human BCL2L1/BCL2L/BCLX (His tag)
recombinant protein