

Mouse Lgmn/Prsc1 (His tag) recombinant protein



Catalog Number: 518354

General Information

Gene Name Synonym

Legumain;Lgmn;Asparaginyl endopeptidase;Protease cysteine 1;Prsc1

Protein Construction

A DNA sequence encoding the Mouse Lgmn/Prsc1 Val18-Tyr435 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris,150mM NaCl,20% Glycerol,pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Lgmn/Prsc1 (His tag) recombinant protein