

# E.coli rppH/nudH/ygdP/b2830/JW2798 recombinant protein



Catalog Number: 518361

## General Information

### Gene Name Synonym

RNA pyrophosphohydrolase, (Di)nucleoside polyphosphate hydrolase, Ap5A pyrophosphatase, rppH

### Protein Construction

A DNA sequence encoding the E.coli rppH/nudH/ygdP/b2830/JW2798 Met1-Gly176 is expressed.

### Organism

E.coli

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 50mM Tris, 500mM NaCl, 10% glycerol, pH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



E.coli rppH/nudH/ygdP/b2830/JW2798 recombinant protein