

# Human CPA4/CPA3/UNQ694/PRO1339 (His tag) recombinant protein



Catalog Number: 518366

## General Information

### Gene Name Synonym

Carboxypeptidase A4, Carboxypeptidase A3, CPA4, CPA3

### Protein Construction

A DNA sequence encoding the Human CPA4/CPA3/UNQ694/PRO1339 Gly17-Tyr421 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human CPA4/CPA3/UNQ694/PRO1339 (His tag)  
recombinant protein