

Mouse Vsig8 (His tag) recombinant protein



Catalog Number: 518368

General Information

Gene Name Synonym

V-set and immunoglobulin domain-containing protein 8; Vsig8

Protein Construction

A DNA sequence encoding the Mouse Vsig8 Val22-Gly262 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, pH7.4, 10% Glycerol.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Vsig8 (His tag) recombinant protein