

Human IMPDH2/IMPD2 (His tag) recombinant protein



Catalog Number: 518401

General Information

Gene Name Synonym

Inosine-5'-monophosphate dehydrogenase 2, IMPDH2, IMP dehydrogenase 2, IMPD 2, IMPDH 2, IMPDH-II, IMPD2.

Protein Construction

A DNA sequence encoding the Human IMPDH2/IMPD2 Met1-Phe514 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris 150mM NaCl pH 8.0 .

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human IMPDH2/IMPD2 (His tag) recombinant protein