

# Human PITPNA/PITPN (His tag) recombinant protein



Catalog Number: 518425

## General Information

### Gene Name Synonym

Phosphatidylinositol Transfer Protein Alpha Isoform, PI-TP-Alpha, PtdIns Transfer Protein Alpha, PtdInsTP Alpha, PITPNA, PITPN

### Protein Construction

A DNA sequence encoding the Human PITPNA/PITPN Met1-Asp270 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM EDTA, 1mM DTT, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PITPNA/PITPN (His tag) recombinant protein